

DEBUG AQUATIC CHAIR



DEBUG AQUATIC CHAIR



The Debug Aquatic Chair is great for pools, waterparks and showers and can even be used as a manual street wheelchair. As a pool wheelchair it can assist in bringing commercial pools with ramps and waterparks into compliance with the new ADA Swimming Pool Guidelines. The interchangeable front axle and sleeves, along with the rear pool suspension, means it only takes a few moments to convert the Debug Aquatic Chair from one style to another.

LONG LASTING #316L STAINLESS STEEL FRAME / REMOVE-ABLE ARTICULATING TELESCOPIC LEG REST / FITS THROUGH A STANDARD ADA DOORWAY / TUBELESS SOLID CORE TIRES / POLYCARBONATE MAG WHEELS ROLL ACROSS GRAVEL / ALL STAINLESS STEEL AND RUST PROOF AXLE AND FASTENERS / DUAL STITCHED SLING WITH TENERA 130 THREAD / PHIFERTEX PLUS SLING MATERIAL RESISTS STRETCHING, MILDEW AND BACTERIA, COLD CRACK AND ULTRAVIOLET DAMAGE AND IS MACHINE WASHABLE.

SPECIFICATIONS

SEAT SIZES	ENGLISH (INCH)	METRIC (CM)
Widths	14, 18, 22	36, 46, 56
Average depth	17	43
Seat-to-floor height	17 - 21	43 - 53
Backrest height	18	46
Telescope footrest <i>knee to heel</i>	15 - 22	38 - 56
Aquatic chair width	21	53

WEIGHTS

Capacity	250 - 400 lbs	114 - 182 kg
Chair weight	35 lbs	16 kg

Manufactured in the USA,
Debug has a 10 year warranty on stainless steel frames and a 30 day money back guarantee.

© Copyright Deming Designs Inc., 2001 - 2016. All rights reserved.

ORDER ONLINE WWW.DEBUGMOBILITYPRODUCTS.COM