

INGRID




In 1998, Popab consulted with leading Swedish experts from two separate industries: design and steel. The aim was to develop the ultimate stick and crutch holder. The result was INGRID, a unique combination of elegance and function. The INGRID holder is appealing to the eye, but its strength lies just as much in its usefulness and functionality as a stick and crutch holder. This is probably why INGRID accessorizes hundreds of thousands of sticks and crutches around the world.

HOOK



PROP



SIMPLE SNAP-ON FITTING



FUNCTION

With the INGRID stick and crutch holder, it is easy to hang your stick or crutch when not in use. Since INGRID has both a hook and a prop, it is possible to hang your stick or crutch on everything from chairs and table ledges to shopping carts and countertops.

DESIGN

The INGRID stick and crutch holder is fashioned in cooperation with leading Swedish ergonomic designers to guarantee user performance and satisfaction. The holder is easy to handle even for users who have reduced muscular strength. Thanks to its exceptional design, INGRID only weighs 26 grams (less than 1 oz). INGRID is made of high quality thermoplastic rubber (Santoprene)

which, together with the stainless steel spring developed especially for INGRID, makes it highly durable. INGRID has a simple snap-on fitting and is suitable for sticks and crutches with various dimensions: 19–25 mm ($\frac{3}{4}$ – 1 inch). A smaller version is available called MINI-INGRID stick holder that fits narrower sticks: 12–16 mm ($\frac{1}{2}$ – $\frac{5}{8}$ inch).

POPAB

Popab's concept has always been to provide walking aid accessories of the highest quality. In addition to the two holders, Popab manufactures over 10 different kinds of cane and crutch tips (ferrules). Popab's products are made in Stockholm, Sweden, using an advanced injection molding technique. The products are pa-

tented and conform to CE directives for Medical Devices.

The product range is currently available on four continents in over 25 countries including the U.S., the United Kingdom, France, Italy, Japan, and Australia. Sales channels include pharmacies, retailers, specialty shops, public healthcare, and manufacturers of walking sticks and crutches.

CONTACT

Postal address: Västermovägen 78
SE-125 40 Älvsjö, Sweden
Phone: +46 (0) 8991220
Fax: +46 (0) 8990515
E-mail: info@popab.se
www.popab.se