

Made in U.S.A.

GREAT STUFF!
Scooters
Wheelchairs
For & Walkers



Diestco Manufacturing Corporation • P.O. Box 6504 • Chico, CA 95927 • 1-800-795-2392

Small Oxygen Tank Holder for manual wheelchair

Part No. B6321

Now with 300 yard visible reflective safety tape!!



A very safe way to carry your oxygen cylinder
on the back of your manual wheelchair.

Products Proudly Made in the U.S.A.

 **DIESTCO Mfg. Corp. • 800-795-2392 • www.diestco.com**


Small O2 Tank Holder For Manual Wheelchairs

Attaches to the back of your manual wheelchair to carry your oxygen cylinder. Made from high quality, water resistant polyester.

DIRECTIONS FOR ATTACHING YOUR SMALL O2 TANK HOLDER TO YOUR WHEELCHAIR.



Slip top loop over wheelchair handle and pull tight. Repeat with second loop.



After bag is secure; the oxygen cylinder can be placed inside. Pull the locking cord tight to ensure that tank is secure.



Pictured above is a centered O2 Tank Holder. Straps can be adjusted to position holder to the left or the right to accommodate your preference.

OTHER EQUIPMENT HOLDERS AVAILABLE FOR YOUR SCOOTER OR POWERCHAIR:

Crutch Holder: Fits forearm and underarm style crutches. Goes behind your scooter or powerchair seat so you can travel with no interferences.

Cane Holder: Allows cane to be stored on back of vehicle when not in use.

Walker Holder: Attaches to your scooter seatback to carry your basic style walker wherever you go!