

## Cane/Sling Seat

HCPCS: E0105



FIGURE A



FIGURE B

- Provides a comfortable seat to rest on when open, and a sturdy quad cane when closed (Figure A)
- Manufactured with sturdy, extruded aluminum tubing
- Four sturdy legs with vinyl contoured tips
- Nylon sling bench, 15" x 5.75"
- Foam grip handles
- 34" tall when closed
- Seat to floor height: 20.5"
- Weighs 2.5 lbs
- 250 lb Weight Capacity

### PRODUCTS

### Description

RTL10360	Cane Sling Seat,Folding,1 EA
----------	------------------------------

### SPECIFICATIONS

### Description

Warranty	Limited Lifetime on frame, 6 Month Warranty on Sling
Weight Capacity	250 lbs
Seat To Floor Height	20.5"
Dimensions	15" (W) x 5.75" (H)