

Nitro Sprint Rollator

HCPCS:



The sleek Nitro Sprint Rollator is designed to provide comfort, convenience, and control at every turn. The rollator sports a unique braking system that reduces tire wear and has an integrated slow-down brake for safer maneuverability over diverse terrain. When not in use, the rollator can be folded side to side with one hand and locked in place for easy transport. Other features include multiple storage bag mounting points and a wider backrest with a window for improved comfort and visibility.

- Key Features & Benefits
- Contoured backstrap with viewing window for improved comfort and visibility.
- The patented frame design provides added strength and durability with a 350lb weight capacity.
- Integrated slow-down brake controls speed when navigating diverse terrain.
- The enhanced braking system reduces tire wear and allows for easy replacement of brake pads and cables.
- Improved tire tread for better stopping ability.
- Built-in accessory mount easily attaches a cup, cane or smartphone holder.
- Built-in brake cable for improved durability and safety
- Large 10" front wheels provide optimal stability and rolling comfort.
- Folds side to side with one hand and remains locked for easy transport and storage.
- A padded seat provides added comfort while sitting.
- Storage bag mounts in multiple locations for easy access.

PRODUCTS	Description
102662BK	Nitro Sprint Rollator,Black1EA
102662BK-H	Nitro Sprint Rollator,Hemi,Black,1EA
102662BK-T	Nitro Sprint Rollator,Tall,Black,1/ea
102662BL	Nitro Sprint Rollator,Blue1EA
102662BL-H	Nitro Sprint Rollator,Hemi,Blue,1EA
102662BL-T	Nitro Sprint Rollator,Tall,Blue,1/ea
102662RD	Nitro Sprint Rollator,Red,1EA
102662RD-H	Nitro Sprint Rollator,Hemi,Red,1EA
102662RD-T	Nitro Sprint Rollator,Tall,Red,1/ea

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	19 lbs
Handle Height	Hemi: 29-33 Standard: 30-36.5 Tall: 32-39.5
Weight Capacity	350 Lbs
FireRating	Back and Seat fabric TB117-2013 section 1
Seat To Floor Height	Standard: 21" Tall: 23 Hemi: 18

SPECIFICATIONS**Description**

Dimensions Assembled	24" (W) X 27.5" (L) X (30" - 36.5") H
Width inside back legs	19"

ACCESSORIES**Description**

102662-CUP	Nitro Universal Cup Holder, 1EA, EA
102662-CANE	Nitro Cane Holder, 1BX, EA
102662-CELL	Nitro Universal Cell Phone Holder, 1EA, EA
102662-TRAY	Nitro Sprint Rollator Tray, 1EA, EA