



Purpose

Our Basic Foam Cushion features a multi-layer design for all-day comfort and support. The one-inch memory foam top molds to your body for customized comfort, while the supportive foam base ensures stability. The cushion is covered in a water-resistant, 2-way stretch Low Shear II fabric, offering skin protection by minimizing friction and easing cleaning for long-lasting use.

Benefits include:

- **Pressure Distribution and Support:** Dual foam layers balance soft comfort with firm support for optimal stability while evenly distributing weight and minimizing pressure.
- **Enhanced Stability:** Prevents shifting to help maintain proper posture and balance.
- **Skin Protection:** Low-shear cover reduces friction and prevents skin irritation.
- **Lightweight and Easy Handling:** Vinyl-free design makes it easier to lift and reposition.

Easy to Use

Place the non-skid fabric bottom on the seating surface and secure with the straps and side-release buckle for a secure fit.

Easy Care

The cover can be wiped clean with a damp cloth and a mild cleanser or disinfecting wipes.

	2" Foam Base	3" Foam Base	4" Foam Base
Size	ReOrder Number	ReOrder Number	ReOrder Number
16"W x 16"D	752001	752002	752003
18"W x 16"D	752004	752005	752006
20"W x 16"D	752007	752008	752009
18"W x 18"D	752010	752011	752012
20"W x 18"D	752013	752014	752015
22"W x 18"D	752016	752017	752018
24"W x 18"D	752019	752020	752021