Popab INGRID



## **INGRID** Ferrules

Ingrid high quality ferrules are not only extremely durable but do also provide for an excellent grip both indoors and outdoors.

Thanks to the tread pattern at the bottom of the ferrule not wearing off as easily as a conventional ferrule may, safety is enhanced.

Three distinct surfaces: Classic, Ice Spike and Indoor







Keep retracted for use during the summer months or while indoor so not to damage floors



The elastic neck adjusts to the dimension of the cane or crutch shaft



Available in 15 combinations of style and size to fit virtually any cane or crutch

## Three distinct surfaces

- Ingrid Classic: The tire track inspired treads at point of contact provides an excellent grip on most surfaces.
- Ingrid Retractable Spike: The retractable spike provides a firm grip on snow, ice and most other outdoor surfaces.
- Ingrid Indoor: For users needing additional support from their cane or crutch. For indoor use.

## Firm or elastic neck

All ferrules have either a firm neck in line with traditional ferrules or the patented elastic neck which, thanks to its design, adjusts to the dimension of the cane or crutch shaft and makes it easy for the user to attach and remove the ferrules as needed.

## Available in 15 combinations of style and size

	Elastic neck				Firm neck	
	Fitness walking poles	Thin canes	All custom canes, crutches	Very thick canes, crutch- es	Thin custom canes, crutches	Many custom canes, crutch- es
	8-12 mm	13-17 mm	16-22 mm	21-27 mm	16-18 mm	19-21 mm
Classic						
Ice spike						
Indoor				- August - A		

The Ingrid products are made in Sweden www.popab.se

