

GREAT STUFF!

Scooters
Wheelchairs
& Walkers

Made in U.S.A.



FROM

DIESTCO
MANUFACTURING
CORPORATION

Heavy Duty Scooter Cover Regular

Part # V1111



600x300 polyester with PVC rubber backing provides the highest level of protection from the elements for your scooter. Works great on lifts or for just storing outside.

Made in the U.S.A.

Heavy Duty Scooter Cover - Regular

V1111

The cover will provide protection for your vehicle from the weather and dust. It is lightweight and easy to get on and off.

Slip the cover over your vehicle and pull the elastic down under your chassis. Pull the 1" Webbing Belt so that it is tight around the bottom. For extra security, when transporting your vehicle on the back of your car or van, tie down the bottom by hooking bungee-type cords (not included) through the bottom belt loops.

- This is a Heavy Duty cover, you will have a 1" webbing belt around the base to tighten instead of a cord lock.
- Depending on your vehicle, cover may not fit over wheels.
- The cover is not intended for highway use without using some securing device (i.e. rope or bungee cords) over the top of the cover. This will reduce the amount of excess room under the cover and prevent the cover from billowing up at high speeds. If your cover is flapping excessively at high speeds, the life of the cover will be greatly reduced.



DIESTCO Mfg. Corp.
P.O. Box 6504
Chico, CA 95927
800-795-2392
www.diestco.com

We strongly encourage you to call us at **800-795-2392** if you have any problems, questions, or comments. We are always developing new products for wheelchairs, scooters, crutch, and walker users...so let us know if there's something special you need.